

92506 Ambulatory Blood Pressure Report Management System



Product Specifications

Interfaces to the 90207 and 90217 Ambulatory Blood Pressure Monitors

Windows-based, mouse-driven software for ease of use

Simple Windows Wizards interface for initializing and uploading information from ABP monitors

Facilitates selection of measurement intervals, inflation frequencies, display of graphic and tabular data, data analysis, and data transfer

Choice of multiple, user-configurable report formats and commonly used data formats to simplify redundant entries

Supports importing of data from 90121 analysis software

Supports exporting of data, such as raw data readings, statistical summaries, and patient information, to ASCII data files

Supports faxing of reports (requires Microsoft Fax or third-party fax utility)

Supports e-mail transfer of reports as PDF documents (requires Adobe Acrobat Writer application and Microsoft Exchange or Outlook Express as e-mail transport mechanism)

Supports database for storing information

Supports multiple users on a single PC system

Supports direct application access to operator manuals


Supports compliance with FDA 21 CFR Part 11 — electronic records, electronic signatures, and the Health Insurance Portability and Accountability Act (HIPAA) regulations Privacy Rule

Supports user logon and password for controlling access to information

Standard Components

Software CD (P/N 068-0230-xx)

Customer Installation Instructions (P/N 075-0422-xx)



92506 Ambulatory Blood Pressure Report Management System



Computer Support

Hardware Minimum Configuration

Pentium CPU 200 MHz or higher

64 MB of memory

CD-ROM drive

Available serial port or USB port

Two- or three-button mouse

SVGA monitor (1024 × 768 resolution or higher)

Required hard disk space

ABP Report Management System Client Application
12 MB

Microsoft Desktop Engine (MSDE)
75 MB

Electronic Manuals (optional)
2 MB

Free hard drive space for Windows virtual memory
and ABP Client Application Database expansion
(additional 50 KB for each ABP test) 250 MB

Software Minimum Configuration

Microsoft Windows 98, 2000, or XP

Printer Support

Any printer supported by Windows with a minimum resolution of 150 dpi

Local or network-based

Required Accessories

Direct serial cable for the 90207/90217 monitors
(P/N 012-0097-02)

Regulatory Approvals



CE marked in accordance with the Medical Device Directive, 93/42/EEC

Other brands and product names are trademarks of their respective owners.

www.spacelabs.com

061-1283-00 Rev. B 09/2004